

1.4QT MAIL-IN TAKE BACK INSTRUCTIONS

1. Close the container by snapping on the permanent locking lid. After the lid is securely closed, apply the security tape strip across the top of the lid.
2. Place container inside the zipper bag. Remove as much air as possible.
3. Place both the container and bag assembly within the ship-back box. Close the inner box flaps first; the shorter flaps needs to close first.



Step 1



Step 2



Step 3

4. Partially close the outer box.
5. Close the two opposing adhesive backed flaps, peel adhesive liner, then firmly shut the non-adhesive backed flaps - sealing the outer box closed.



Step 4



Step 5

6. Bring your filled and sealed ship-back box to your local UPS® office to ship.

Only for disposal of materials free of pathogens and transferable diseases. If waste was exposed, it must be disinfected prior to disposal. Drain all free liquids prior to disposal.
DO NOT OVERFILL.

Track your product online at www.americanrxgroup.com/tracking/
Need help? Call 877-350-0693 or visit www.americanrxgroup.com